

Fully documented source code for the ReporterIonAgent

```
1  package peptizer.agents.custom;
2
3  import be.proteomics.mascotdatfile.util.interfaces.Spectrum;
4  import be.proteomics.mascotdatfile.util.mascot.Peak;
5  import be.proteomics.mat.MatConfig;
6  import be.proteomics.mat.interfaces.Agent;
7  import be.proteomics.mat.util.AgentReport;
8  import be.proteomics.mat.util.PeptideIdentification;
9  import be.proteomics.mat.util.enumerator.AgentVote;
10 import be.proteomics.mat.util.fileio.MatLogger;
11
12 import java.math.BigDecimal;
13 import java.util.Properties;
14 /**
15  * Created by IntelliJ IDEA.
16  * User: kenny
17  * Date: 18-jun-2008
18  * Time: 15:30:28
19  */
20
21 /**
22  * Class description:
23  * -----
24  * This class was developed to inspect for deviating reporter ion intensities.
25  */
26 public class ReporterIonAgent extends Agent {
27
28     /**
29      * PARAMETERS
30      * -----
31      * String identifiers for the parameters in the agent.xml configuration file.
32      */
33
34     // Mass over charge parameter for the first reporter ion.
35     public static final String MASS_1 = "reporter_mz_1";
36
37     // Mass over charge parameter for the second reporter ion.
38     public static final String MASS_2 = "reporter_mz_2";
39
40     // Fold ratio parameter between the two reporter ions you consider as deviating.
41     public static final String RATIO = "ratio";
42
43     // Error tolerance for matching the expected reporter
44     // ion mass over charge values to a fragmentation in the MS/MS spectrum.
45     public static final String ERROR = "error";
46
47     /**
48      * CONSTRUCTOR
49      * -----
50      * Construct a new instance of the ReporterIonAgent.
51      */
52     public ReporterIonAgent() {
53         // 1. Calls the constructor of the Agent superclass.
54         super();
55         // 2. Gets the properties for this Agent.
56         // A singleton MatConfig object reads the agent.xml configuration file upon starting Peptizer.
57         // The Properties of each Agent can be then be retrieved by the unique identifier of an Agent.
58         Properties prop = MatConfig.getInstance().getAgentProperties(this.getUniqueID());
59
60         try {
```

Fully documented source code for the ReporterIonAgent

```
61      // 2a. Sets common properties shared among all Agents.
62      super.setName(prop.getProperty("name"));
63      super.setActive(Boolean.valueOf(prop.getProperty("active")));
64      super.setVeto(Boolean.valueOf(prop.getProperty("veto")));
65
66      // 2b. Sets specific properties for this ReporterIonAgent.
67      // Masses of the two reporter ions, the threshold ratio and the mass error tolerance.
68      this.iProperties.put(MASS_1, prop.getProperty(MASS_1));
69      this.iProperties.put(MASS_2, prop.getProperty(MASS_2));
70      this.iProperties.put(RATIO, prop.getProperty(RATIO));
71      this.iProperties.put(ERROR, prop.getProperty(ERROR));
72  } catch (NullPointerException npe) {
73      // Note that an exception is thrown when one of the parameters is missing!
74      MatLogger.logExceptionalGUIMessage("Missing Parameter!!",
75          "Parameters " + MASS_1 + ", " + MASS_2 + ", " + RATIO + "and " + ERROR +
76          " are required! for Agent \"" + this.getName() + "\" !!\nExiting..");
77      System.exit(0);
78  }
79  }
80
81  /**
82   * INSPECTION
83   * -----
84   * The inspection is the core of an Agent since this logic leads to the Agent's vote.
85
86   * This method returns an array of AgentVote objects, reflecting this Agent's idea
87   * whether to select or not to select the peptide hypothesis.
88
89   * All Agent Implementations must also create and store AgentReport for each peptide hypothesis.
90   *
91   * @param aPeptideIdentification PeptideIdentification that has to be inspected.
92   * @return AgentVote[] as a vote upon inspection for each the confident peptide hypotheses.
93   * AgentVotes[0] gives the inspection result on PeptideHit 1
94   * AgentVotes[1] gives the inspection result on PeptideHit 2
95   * AgentVotes[n] gives the inspection result on PeptideHit n+1
96
97   * Where the different AgentVotes can be:
98
99   * a vote approving the selection of the peptide hypothesis.
100  * a vote indifferent to the selection.
101  * a vote objecting to select the peptide hypothesis.
102  */
103  public AgentVote[] inspect(PeptideIdentification aPeptideIdentification) {
104
105      // A. PREPARING THE VARIABLES
106      // *****
107
108      // 1. The reporter ion masses.
109      double lReporterMass_1 = Double.parseDouble((String) (this.iProperties.get(MASS_1)));
110      double lReporterMass_2 = Double.parseDouble((String) (this.iProperties.get(MASS_2)));
111
112      // 2. The fold ratio threshold.
113      double lRatio = Double.parseDouble((String) (this.iProperties.get(RATIO)));
114
115      // 3. The error tolerance.
116      double lError = Double.parseDouble((String) (this.iProperties.get(ERROR)));
117
118      // 4. Reserves an array with AgentVotes for each confident peptide hypothesis.
119      AgentVote[] lAgentVotes =
120          new AgentVote[aPeptideIdentification.getNumberOfConfidentPeptideHits()];
```

Fully documented source code for the ReporterIonAgent

```
121
122     // Since this inspection is dependent on the MS/MS spectrum,
123     // it will result in the same vote for each peptide hypothesis.
124     // Therefore, a single inspection is reused for each peptide hypothesis.
125
126     // 5. Initiate an AgentReport serving as a report for this inspection.
127     iReport = new AgentReport(getUniqueID());
128
129     // 6. Local variable to store the result shown in the information table.
130     String lResultForTable = "";
131     // 7. Local variable to store the result written in the arff file.
132     String lResultForArff = "";
133
134     // 8. Local variables for matching reporter ion 1 in the MS/MS spectrum.
135     boolean lReporter_1_match = false;
136     double lReporter_1_intensity = 0;
137
138     // 9. Local variables for matching reporter ion 2 in the MS/MS spectrum.
139     boolean lReporter_2_match = false;
140     double lReporter_2_intensity = 0;
141
142     // B. THE ACTUAL INSPECTION
143     //*****
144
145     // 1. Gets the MS/MS spectrum from the PeptideIdentification object
146     // that was given as a parameter to the inspect() method.
147     Spectrum lSpectrum = aPeptideIdentification.getSpectrum();
148
149     // 2. Gets the peaklist from this MS/MS spectrum.
150     Peak[] lPeaks = lSpectrum.getPeakList();
151
152     // 3. Iterates over all peaks through a for loop.
153     for (int i = 0; i < lPeaks.length; i++) {
154         Peak lPeak = lPeaks[i];
155         double lDelta_1 = lPeak.getMZ() - lReporterMass_1;
156         double lDelta_2 = lPeak.getMZ() - lReporterMass_2;
157
158         // 3i) If absolute value of the mass difference of this fragmentation and the expected mass
159         // of reporter ion 1 is less then the defined error tolerance.
160         // Then there is a match!
161         if (Math.abs(lDelta_1) < lError) {
162             lReporter_1_match = true;
163             lReporter_1_intensity = lPeak.getIntensity();
164         }
165
166         // 3ii) Same idea for reporter ion 2.
167         if (Math.abs(lDelta_2) < lError) {
168             lReporter_2_match = true;
169             lReporter_2_intensity = lPeak.getIntensity();
170         }
171
172         // 3iii) For performance reasons: if both peaks were matched then the for loop can be exited.
173         if (lReporter_1_match && lReporter_2_match) {
174             break;
175         }
176     }
177
178     // 4. Checks the intensity ratio.
179     double lUpperBoundary;
180     double lLowerBoundary;
```

Fully documented source code for the ReporterIonAgent

```
181
182 // First an upper and a lower boundary must be defined for the Reporter Ions intensity ratio.
183 // If the experimental ratio between Reporter Ion 1 and Reporter Ion 2 is
184 // more then 1.5 or less then 0.66, then the two samples deviate by a factor of 1.5.
185 //
186 // If the user defined the factor as larger then 1, the upper boundary for an
187 // deviating ratio is given by 1 divided by that factor.
188 // The lower boundary for an deviating ratio is given by 1 multiplied by that factor.
189 // Example:
190 // If (Ratio=1.5)
191 // Then lower boundary = 1.5 and upper boundary = 0.66
192 if (lRatio > 1) {
193     lUpperBoundary = 1 * lRatio;
194     lLowerBoundary = 1 / lRatio;
195 } else {
196     // If the user defined the facotr as smaller then 1, then it is the other way round.
197     lUpperBoundary = 1 / lRatio;
198     lLowerBoundary = 1 * lRatio;
199 }
200
201 // 5. Local variable for the intesity ratio between Reporter Ion 1 and Reporter Ion 2.
202 double lExperimentalRatio;
203 // Local boolean for the upcoming function to store whether the ratio between the Reporter Ions
204 // deviates more then the expected ratio.
205 boolean lDeviatingRatio;
206
207 // 6. Checks the intensities of the reporter ions!
208
209 // 6i) If this condition is true, then one of the reporter ions was not found!
210 // These are not selected as these are probably unlabeled peptides.
211 if (!lReporter_1_match || !lReporter_2_match) {
212     lDeviatingRatio = false;
213     lExperimentalRatio = 0;
214     // 6ii) Else both the reporter ions were found. Lets inspect their ratio.
215 } else {
216     lExperimentalRatio = lReporter_1_intensity / lReporter_2_intensity;
217
218     // The ExperimentalRatio between reporter ion 1 an reporter ion 2
219     // is either less then the lower boundary,
220     // or either more then the upper boundar.
221     // In both cases, the reporter ion intesity is deviating for both samples:
222     if (lLowerBoundary > lExperimentalRatio || lUpperBoundary < lExperimentalRatio) {
223         // A. The Agent inspection resulted in deviating reporter ion intensities as
224         // their ratio was outside one of the boundaries.
225         lDeviatingRatio = true;
226     } else {
227         // B. Else the Agent inspection resulted in non deviating reporter ion intensities as
228         // their ratio was within the lower and upper boundary.
229         lDeviatingRatio = false;
230     }
231 }
232
233 // C. MAKING THE INSPECTION REPORTS AND COMMITTING THE VOTES
234 //*****
235
236 // 1. In all cases, store the experimental ratio between the
237 // two reporter ions as a value to display in the information table.
238
239 // A BigDecimal rounds a double at 2 decimals.
240 BigDecimal lRoundedExperimentalRatio = null;
```

Fully documented source code for the ReporterIonAgent

```
241     lRoundedExperimentalRatio = new BigDecimal(lExperimentalRatio).setScale(2, BigDecimal.ROUND
    _HALF_UP);
242     lResultForTable = lRoundedExperimentalRatio.toString();
243
244     // 2i. Deviating reporter ion intensity ratio, this Agent suggests to select the peptide hypothesis!
245     if (lDeviatingRatio) {
246         lResultForArff = "1";
247         for (int i = 0; i < lAgentVotes.length; i++) {
248             lAgentVotes[i] = AgentVote.POSITIVE_FOR_SELECTION;
249         }
250     // 2ii. Non Deviating reporter ion intensity ratio, this Agent is neutral to select the peptide hypothesis!
251     } else {
252         lResultForArff = "0";
253         for (int i = 0; i < lAgentVotes.length; i++) {
254             lAgentVotes[i] = AgentVote.NEUTRAL_FOR_SELECTION;
255         }
256     }
257
258     // 3. Creates an Agentreport for this inspection.
259     iReport.addReport(AgentReport.RK_RESULT, lAgentVotes[0]);
260     iReport.addReport(AgentReport.RK_TABLEDATA, lResultForTable);
261     iReport.addReport(AgentReport.RK_ARFF, lResultForArff);
262
263     // 4. Stores the report on the PeptideIdentification object.
264     for (int i = 0; i < lAgentVotes.length; i++) {
265         aPeptideIdentification.addAgentReport((i + 1), this.getUniqueID(), iReport);
266     }
267
268     // 5. Returns the AgentVotes in the end of the inspection.
269     return lAgentVotes;
270 }
271
272 /**
273  * Returns a description for the Agent.
274  * Note that html tags are used to stress properties.
275  * Use in tooltips and configuration settings.
276  * <p/>
277  * Fill in an agent description. Report on purpose and a minor on actual implementation.
278  *
279  * @return String description of the ReporterIonAgent.
280  */
281 public String getDescription() {
282     return "<html>Inspects for the aberrant reporter ion intensities." +
283         "<b>Selects when two reporter ions ( " + this.iProperties.get(MASS_1) +
284         " , " + this.iProperties.get(MASS_2) + ") have a more then " + this.iProperties.get(RATIO) +
285         " fold intensity ratio.";
286 }
287 }
```